

PRODUCTS & SERVICES

# cooltech

## FREEZE OUT FAT

TOTAL BODY CONTOURING WITH SIX DIFFERENT APPLICATORS FOR 20-30 PER CENT FAT REDUCTION AFTER ONE TREATMENT

**C**ooltech is a premium fat-freezing technology based on the proven science of cryolipolysis, studied in numerous clinical trials with certified clinical predictability and maximum results.

Distributed in Australia by Cryomed, the European manufactured device is designed for treating localised fat pockets. Not only is Cooltech affordable with low consumables, it also presents a range of unique advantages to produce fast and effective results, while helping secure a strong ROI.

Its six applicators are fashioned to adapt to any area of the body, offering a total body contouring treatment.

Cooltech works by gently drawing fatty skin into applicators, before cooling the area to the desired treatment temperature for maximum clinical effect.

The different hand pieces are adapted for areas of curved contours like the back, waist and sub-buttocks, broad areas such as the flanks, love handles and calves, and for fat folds small and large, from the abdomen and back to the arms, knees and under the chin.

The procedure targets and cools fat cells to temperatures that trigger fat cell apoptosis and elimination without affecting any of the other surrounding tissue.

All the hand pieces work at a temperature between 3°C and -8°C. This temperature range is required to treat fat folds of various thicknesses.

The Cooltech device allows two handpieces to be applied simultaneously, thereby optimising treatment times and saving you and your patients valuable time.

The full set of handpieces allows operators to provide full body treatments to their patients with results comparable to, if not better than, any other registered cryolipolysis device on the market.

### FAST AND EFFECTIVE RESULTS

After one treatment patients can expect to see a 20-30 per cent fat reduction in the treated area. Results can first be seen in eight to 12 weeks.

Typically, patients are satisfied with the results of the first treatment but, if necessary, repeat treatments can be performed at three months, with a maximum of three sessions per area.

Cooltech involves no needles, incisions, downtime or compression garments. It's an easy and reliable treatment that allows the patient to immediately resume normal activities.

### THE APPLICATORS

#### Straight HP

Suitable for areas such as the abdomen, inside of the knees, sides, etc.

#### Curved HP

Indicated for fat located on curved contours. It works on areas like fat folds on the back, the waist, sub-buttocks, etc.

#### Tight HP

Works on areas such as the abdomen, sub-buttocks, back, hips, etc.

#### Double HP

Designed exclusively for treating the abdominal area.

#### Oval HP

Works on sides, love handles, calves, etc.

#### Tiny HP

Works on areas such as arms, sub-scapula, sub-umbilical accumulations and submentum.

### INTERNATIONAL CERTIFICATIONS, PROVEN RESULTS

Now available in more than 40 countries, over one million Cooltech treatments have been performed worldwide.

Numerous clinical studies have demonstrated its effectiveness, along with very high satisfaction levels across a wide variety of patients.

The Cooltech device complies with the European CE 0051 medical certificate, in addition to the ISO 9001 quality certification and ISO 13485 quality standard. In Australia, it is entered in the Australian Register of Therapeutic Goods (ARTG), inclusion no: 276394.

Cooltech is manufactured by Cocoon Medical, a Spanish company with its head office in Barcelona.

The company's highly skilled researchers and innovators have produced equipment that has gained worldwide attention for its innovation and excellence.

Cocoon Medical's clinical unit helps to ensure the effectiveness and safety of all Cooltech devices on the market, and its clinical monitoring and training unit supports both the in vitro and in vivo therapeutic excellence of all devices.

Distributor Cryomed provides a 5-star level of service and ongoing support for its clients, ensuring optimal training and education to help yield the best possible results for the patient and the clinic. **AMP**

For more information, call Cryomed on 1300 346 448 or visit [www.cryomed.com.au](http://www.cryomed.com.au)

### EXCEEDING EXPECTATIONS

Dr Jenny Kimmins, of Southern Cosmetics in Melbourne, says the clinic has been using Cooltech since January and 'patient results overall have exceeded our expectations.'

'We see on average two to three Cooltech patients a week and are so confident of its capabilities that we as good as guarantee our treatments.'

Dr Kimmins says the lower abdomen/stomach is the most highly requested area for treatment, but the clinic is finding that patients are so happy with the outcomes achieved, they come back to have their upper abdomen treated as well.

Other areas that have proved popular are love handles/flanks and the pockets of fat beneath the bra line that concern so many women.

Fleshy upper arms is another area that Dr Kimmins says has reaped surprisingly positive results.

'When a patient has a lot of loose skin as well as fat on their upper arms, surgery is usually the only option – other than disguising them with clothing or just accepting the fact,' she says. 'This is the case when, say, a patient has lost a lot of weight, creating redundant skin.'

'When we assess patients at the initial consultation, whatever area of treatment they are seeking, we will always say if we believe surgery, such as liposuction, is a better option to achieve their desired outcome.'

'But if a patient is insistent they want to try Cooltech anyway, we will consider their wishes.'

'Recently, a 29-year-old woman came to us with a lot of loose skin and residual fat under her arms after a significant weight loss. She decided to proceed with Cooltech and has had a spectacular result.'

Generally speaking, Dr Kimmins says an ideal candidate for Cooltech is a person of average or slightly above average weight with stubborn fat deposits on certain areas of their body.

